Bardex™ I.C.

Anti-infective Foley Catheter with Bard®
Hydrogel and Bacti-Guard®* Silver Alloy
Coating

Instructions for Use

en	English	2
fr	French	5
de	German	9
it	Italian	13
es	Spanish	17
nl	Dutch	21
pt	Portuguese	25
el	Greek	29
da	Danish	33
sv	Swedish	37
fi	Finnish	41
no	Norwegian	45
pl	Polish	48
tr	Turkish	52
	Symbols	55

Bardex™ I.C.

Anti-infective Foley Catheter with Bard® Hydrogel and Bacti-Guard®* Silver Alloy Coating

The clinical benefits of Bardex™ I.C. Latex Foley Catheters are that they aid in clinical management by allowing for continuous drainage and measurement of urine. In addition, Bardex™ I.C. Foley Catheters benefit users by lowering the occurrence of bacterial adhesion to the catheter. The hydrogel coating is designed to reduce friction in catheter insertion.

Foley catheters can be used in patients of all ages. Up to 10 CH/Fr is generally considered pediatric sizing; however, the appropriate size should be determined by the healthcare professional.

Recommended indwelling time not to exceed 28 days.

Indications for Use: The Bardex™ I.C. Latex Foley Catheter is indicated for use in the drainage and/or collection and/or measurement of urine. Generally, drainage is accomplished by inserting the catheter through the urethra and into the bladder. However, drainage is sometimes accomplished by suprapubic or other placement of the catheter, such as nephrostomy tract.

Intended User: This product is intended to be used by qualified Healthcare Professionals.

Intended Use: For urological use only.

Contraindications: No known contraindications.

Warnings:

- On catheter, do not use ointments or lubricants having a petroleum base. They will damage catheter and may cause balloon to burst.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.
- This is a single use device. Do not resterilize any portion of this device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness, or death of the patient.

Users and/or patients within the European Union should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users outside of the European Union should report any serious incident that has occurred in relation to the device to the manufacturer and the regulatory authority of the country in which the user and/or patient is established.

Precautions:

- Do not use device if package is opened or damaged.
- Do not aspirate urine through drainage funnel wall.
- Should balloon rupture occur, care should be taken to assure that all balloon fragments have been removed from the patient.

Cautions:

- **This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.**
- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Storage: Keep away from sunlight.

Instructions for Use:

- Follow your facility protocol for all processes related to catheterisation, re-catheterisation, stabilization and urine collection.
- Follow aseptic procedures outlined by your local hospital or healthcare facility. Thoroughly wash hands, wear sterile gloves, and use aseptic technique whenever handling the catheter or catheter components during insertion.
- Visually inspect the product for any imperfections or surface deterioration prior to use. If package is opened or if any imperfection or surface deterioration is observed, do not use.

Insertion:

1. Wash hands.
2. Use aseptic technique to prepare the patient for catheter insertion.
3. Use aseptic technique to remove the catheter from package. Connect catheter to urine collection container as needed.
4. Lubricate the catheter. If hospital policy permits, it is possible to inject lubricant into the urethra.
5. Catheterize the patient using dominant hand.
6. When catheter tip has entered bladder, urine will flow through the catheter.
 - a. For female anatomy, insert catheter 2 more inches (approximately 5 cm).
 - b. For male anatomy, insert catheter all the way to bifurcation.

7. Inflate the balloon with pre-filled water syringe if included.
 - a. If pre-filled syringe is not included, then use a slip tip/luer lock syringe to fill catheter balloon with sterile water. Do not use needle.
 - b. Do not exceed recommended capacities as stated on the applicable labeling of the catheter.
 - c. The chart below provides additional information on inflation volume and balloon size.
 - d. Inflation Volume (A) incorporates the volume required to pass through lumen and fully inflate the specific Balloon Size (B).

	Inflation Volume (A)	Balloon Size (B)
	10 mL	5 mL

Note: Drawing may not display actual shape of catheter tip or presence of additional lumen(s).

8. Gently pull the catheter until the catheter balloon is snug against the bladder neck.

Removal:

1. Wash hands.
2. Gently insert a luer slip tip syringe into the catheter valve.
 - a. Never use more force than is required to make the syringe “stick” in the valve.
3. Allow the pressure within the balloon to force the plunger back and fill the syringe.
 - a. If you notice slow or no deflation, re-seat the syringe gently.
4. Use only gentle aspiration to encourage deflation, if needed. Vigorous aspiration may collapse the inflation lumen, preventing balloon deflation.
5. If balloon does not deflate, contact trained professional for assistance, as directed by hospital protocol.
6. After balloon is fully deflated, remove catheter, and dispose of in accordance with hospital policy.

For U.S. patents that may apply, see [bd.com/patents](https://www.bd.com/patents).

The Becton, Dickinson and Company English symbols glossary is located at [bd.com/en-us/symbols-glossary](https://www.bd.com/en-us/symbols-glossary).

BD, Bard and Bardex are trademarks of Becton, Dickinson and Company or its affiliates.
© 2024 BD. All rights reserved.

*Bacti-Guard® silver alloy coating is licensed from Bactiguard AB.
Bacti-Guard is a registered trademark of Bacti-Guard AB.

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C. R. Bard, Inc.
Becton Dickinson Ireland Ltd.
Becton, Dickinson U.K. LTD