BD[®] Lubri-Sil[™] I.C.

Silicone Foley Catheter with
Bard[®] Hydrogel and Bactiguard^{®*}
Noble Metal Alloy Coating

en	English	2
fr	French/Français	5
de	German/Deutsch	8
it	Italian/Italiano	11
es	Spanish/Español	14
nl	Dutch/Nederlands	17
da	Danish/Dansk	20
sv	Swedish/Svenska	23
fi	Finnish/Suomi	26
no	Norwegian/Norsk	29
pl	Polish/Polski	32
ro	Romanian/Română	35
lt	Lithuanian/Lietuvių k.	38
	Symbols	41

BD[®] Lubri-Sil[™] I.C.

Silicone Foley Catheter with Bard[®] Hydrogel and Bactiguard^{®*} Noble Metal Alloy Coating

The clinical benefits of Lubri-Sil[™] I.C. Silicone Foley Catheters are that they aid in clinical management by allowing for continuous drainage and measurement of urine. In addition, Lubri-Sil[™] I.C. Silicone Foley Catheters benefit users by lowering the occurrence of bacterial adhesion to the catheter. The hydrogel coating is designed to reduce friction in catheter insertion. The Bactiguard[®] noble metal alloy coating is composed of Gold, Silver and Palladium.

Foley catheters can be used in patients of all ages. Up to 10 CH/Fr is generally considered pediatric sizing; however, the appropriate size should be determined by the healthcare professional.

Recommended indwelling time not to exceed 28 days.

Indications for Use: The Lubri-Sil[™] I.C. Silicone Foley Catheter is indicated for use in the drainage and/or collection and/or measurement of urine. Generally, drainage is accomplished by inserting the catheter through the urethra and into the bladder. However, drainage is sometimes accomplished by suprapubic or other placement of the catheter, such as nephrostomy tract.

Intended User: This product is intended to be used by qualified Healthcare Professionals.

Intended Use: For urological use only.

Contraindications: No known contraindications.

Warnings:

- On catheter, do not use ointments or lubricants having a petroleum base. They will damage catheter and may cause balloon to burst.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.
- This is a single use device. Do not resterilize any portion of this device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness, or death of the patient.

Users and/or patients within the European Union should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Users outside of the European Union should report any serious incident that has occurred in relation to the device to the manufacturer and the regulatory authority of the country in which the user and/or patient is established.

Precautions:

- Do not use device if package is opened or damaged.
- Do not aspirate urine through drainage funnel wall.
- Should balloon rupture occur, care should be taken to assure that all balloon fragments have been removed from the patient.

Cautions:

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Instructions for Use:

- Follow your facility protocol for all processes related to catheterisation, re-catheterisation, stabilization and urine collection.
- Follow aseptic procedures outlined by your local hospital or healthcare facility. Thoroughly wash hands, wear sterile gloves, and use aseptic technique whenever handling the catheter or catheter components during insertion.
- Visually inspect the product for any imperfections or surface deterioration prior to use. If package is opened or if any imperfection or surface deterioration is observed, do not use.

Insertion:

1. Wash hands.
2. Use aseptic technique to prepare the patient for catheter insertion.
3. Use aseptic technique to remove the catheter from package. Connect catheter to urine collection container as needed.
4. Lubricate the catheter. If hospital policy permits, it is possible to inject lubricant into the urethra.
5. Catheterize the patient using dominant hand.
6. When catheter tip has entered bladder, urine will flow through the catheter.
 - a. For female anatomy, insert catheter 2 more inches (approximately 5 cm).
 - b. For male anatomy, insert catheter all the way to bifurcation.
7. Inflate the balloon with pre-filled water syringe if included.
 - a. If pre-filled syringe is not included, then use a slip tip/luer lock syringe to fill catheter balloon with sterile water. Do not use needle.
 - b. Do not exceed recommended capacities as stated on the applicable labeling of the catheter.
 - c. The chart below provides additional information on inflation volume and balloon size.
 - d. Inflation Volume (A) incorporates the volume required to pass through lumen and fully inflate the specific Balloon Size (B).

	Inflation Volume (A)	Balloon Size (B)
	10 cc/mL	5 cc/mL
	35 cc/mL	30 cc/mL

Note: Drawing may not display actual shape of catheter tip or presence of additional lumen(s).

8. Gently pull the catheter until the catheter balloon is snug against the bladder neck.

Removal:

1. Wash hands.
2. Gently insert a luer slip tip syringe into the catheter valve.
 - a. Never use more force than is required to make the syringe “stick” in the valve.
3. Allow the pressure within the balloon to force the plunger back and fill the syringe.
 - a. If you notice slow or no deflation, re-seat the syringe gently.
4. Use only gentle aspiration to encourage deflation, if needed. Vigorous aspiration may collapse the inflation lumen, preventing balloon deflation.
5. If balloon does not deflate, contact trained professional for assistance, as directed by hospital protocol.
6. After balloon is fully deflated, remove catheter, and dispose of in accordance with hospital policy.

For U.S. patents that may apply, see **[bd.com/patents](https://www.bd.com/patents)**

The Becton, Dickinson and Company English symbols glossary is located at **[bd.com/en-us/symbols-glossary](https://www.bd.com/en-us/symbols-glossary)**

BD, Bard and Lubri-Sil are trademarks of Becton, Dickinson and Company or its affiliates.
© 2025 BD. All rights reserved.

*Bactiguard® noble metal alloy coating is licensed from Bactiguard AB. Bactiguard is a registered trademark of Bactiguard AB.

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**Consult instructions for use or
consult electronic instructions
for use**

Consulter les instructions d'utilisation
ou les instructions d'utilisation
électroniques • Gebrauchsanweisung
oder elektronische
Gebrauchsanweisung beachten •
Consultare le istruzioni per l'uso in
formato cartaceo o elettronico •
Consulte las instrucciones de
uso o las instrucciones de uso
electrónicas • Raadpleeg de
gebruiksaanwijzing of raadpleeg de
elektronische gebruiksaanwijzing •
Se brugervejledningen, eller se den
elektroniske brugervejledning •
Se bruksanvisningen eller den
elektroniska bruksanvisningen •
Lue käyttöohjeet tai sähköiset
käyttöohjeet • Se bruksanvisningen
eller den elektroniske
bruksanvisningen • Sprawdzić
w instrukcji obsługi lub sprawdzić
w elektronicznej instrukcji obsługi •
Consultați instrucțiunile de utilizare
sau instrucțiunile de utilizare în
format electronic • Susipažinkite
su naudojimo instrukcija arba
elektronine naudojimo instrukcija

**Sterilized using Ethylene Oxide**

Stérilisé à l'oxyde d'éthylène • Mit
Ethylenoxid sterilisiert • Sterilizzato
a ossido di etilene • Esterilizado
mediante óxido de etileno •
Gesteriliseerd met ethyleenoxide •
Steriliseret med ethylenoxid •
Steriliserad med etylenoxid •
Steriloitu etyleenioksidilla •
Sterilisert ved bruk av etylenoksid •
Wysterylizowano za pomocą tlenku
etylenu • Sterilizat cu etilen-oxid •
Sterilizuota etileno oksidu

**Single sterile barrier system**

Système de barrière stérile simple • Einfaches Sterilbarrieresystem • Sistema a barriera sterile singola • Sistema de barrera estéril único • Enkelvoudig steriel barriersysteem • Enkelt, sterilt barriersystem • Enkelt sterilbarriärsystem • Yksinkertainen steriili suojajärjestelmä • Enkelt sterilt barriersystem • System pojedynczej bariery sterylnej • Sistem cu o singură barieră sterilă • Vienkartinė sterili barjerinė sistema

**Do not use if package is damaged and consult instructions for use**

Ne pas utiliser si l'emballage est endommagé et consulter les instructions d'utilisation • Inhalt bei beschädigter Packung nicht verwenden und Gebrauchsanweisung beachten • Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso • No utilizar si el envase está dañado y consultar las instrucciones de uso • Niet gebruiken indien de verpakking beschadigd is en raadpleeg de gebruiksaanwijzing • Må ikke anvendes, hvis emballagen er beskadiget. Se desuden brugervejledningen • Använd inte om förpackningen är skadad och se bruksanvisningen • Ei saa käyttää, jos pakkaus on vaurioitunut, ja lue käyttöohjeet • Ikke bruk dersom pakningen er skadet og se bruksanvisningen • Nie stosować, jeśli opakowanie zostało uszkodzone i sprawdzić w instrukcji obsługi • Nu utilizați dacă ambalajul este deteriorat și consultați instrucțiunile de utilizare • Nenaudokite, jei pakuotė pažeista, ir susipažinkite su naudojimo instrukcija

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Authorized Representative in the European Community

Représentant agréé dans la Communauté européenne • Bevollmächtigter in der Europäischen Gemeinschaft • Rappresentante autorizzato per la Comunità Europea • Representante autorizado en la Comunidad Europea • Gevolmachtigd vertegenwoordiger binnen de Europese Gemeenschap • Autoriseret repræsentant i EU • Auktoriserad representant i Europeiska unionen • Valtuutettu edustaja Euroopan yhteisössä • Autorisert representant i Europa • Autoryzowany przedstawiciel na terenie Wspólnoty Europejskiej • Reprezentant autorizat în Comunitatea Europeană • Įgaliotasis atstovas Europos Bendrijoje

Authorized Representative in Switzerland

Représentant agréé en Suisse • Bevollmächtigter in der Schweiz • Rappresentante autorizzato per la Svizzera • Representante autorizado en Suiza • Gevolmachtigd vertegenwoordiger in Zwitserland • Autoriseret repræsentant i Schweiz • Auktoriserad representant i Schweiz • Valtuutettu edustaja Sveitsissä • Autorisert representant i Sveits • Autoryzowany przedstawiciel w Szwajcarii • Reprezentant autorizat pentru Elveția • Įgaliotasis atstovas Šveicarijoje

Authorized Representative in the United Kingdom

Représentant agréé pour le Royaume-Uni • Bevollmächtigter im Vereinigten Königreich • Rappresentante autorizzato per il Regno Unito • Representante autorizado en el Reino Unido • Gemachtigd vertegenwoordiger in het Verenigd Koninkrijk • Autoriseret repræsentant i Storbritannien • Auktoriserad representant i Storbritannien • Valtuutettu edustaja Yhdistyneessä kuningaskunnassa • Autorisert representant i Storbritannia • Autoryzowany przedstawiciel w Wielkiej Brytanii • Reprezentant autorizat pentru Regatul Unit • Įgaliotasis atstovas Jungtinėje Karalystėje

BARD

C. R. Bard, Inc.
8195 Industrial Boulevard
Covington, GA 30014 USA
+1.844.823.5433
bd.com

Becton Dickinson Ireland Ltd.
Donore Road, Drogheda
Co. Louth, A92 YW26
Ireland

Becton Dickinson
Distribution Center NV
Laagstraat 57, 9140 Temse
Belgium

BD Switzerland Sàrl
Route de Crassier 17
Business Park Terre-Bonne
Bâtiment A4, 1262 Eysins, Switzerland

Becton, Dickinson U.K. LTD
1030 Eskdale Road
Winnersh Triangle
Wokingham, Berkshire
England, RG41 5TS